


Dextran Methacryloyl (DexMA)

Product component

Item	Character	Package Size	Remark
A: DexMA	White spongy	1 g/bottle	Keep in dark
B: Photoinitiator LAP	White powder	0.05 g/bottle	

This instruction applies to EFL-DexMA



DexMA molecular structure

Product introduction

Dextran methacryloyl (DexMA) is a double bond modified dextran, which can be crosslinked and solidified into a gel by UV and visible light under the action of photoinitiator. Due to the excellent water solubility, good biosafety and anti-nonspecific protein adsorption of DexMA, DexMA-based material systems have been widely used in biomedical fields, such as reducing thrombosis in blood vessels, reducing blood viscosity and drug delivery.

Applications

Cell culture, biological 3D printing, tissue engineering, etc.

Storage

Dry kit: room temperature, 3 months; 4°C, 12 months; -20°C, 18 months. **Sterile solution:** 4°C (in dark), 7 days; -20°C (in dark), 6 months. **Please note that repeated freezing and thawing of the solution will affect the performance of the product, so it is best to prepare it when using it.**

Period of validity

The date of manufacture is shown in the package.

Scan the QR code on the right for more information



WeChat Tencent

Solution preparation

1. Prepare 0.25% (w/v) standard solution of initiator

- (1) Add 20ml PBS into the brown bottle containing initiator LAP (containing 0.05g LAP);
- (2) Heat and dissolve the solution in a water bath at 40-50°C for 15 minutes, shaking several times.

The LAP standard solution can be stored for 12 months at 4°C in dark.

2. Prepare DexMA solution (5-15% (w/v) is recommended)

- (1) Take the required mass of DexMA into the centrifugal tube;
- (2) Add the initiator standard solution into the centrifuge tube;
- (3) Dissolve at room temperature for 30 minutes, oscillating several times during the period;
- (4) Sterilize the DexMA solution immediately with a 0.22μm sterile needle filter.

Suggestions for 2D cell culture

- Injecting DexMA solution into the well plate immediately;
(96-well plate: 50-100 μL/ well, 48-well plate: 100-300 μL/ well, 24-well plate: 300-500 μL/ well);
- Irradiate the wells with 405 nm light for 10-30 seconds to gelate, the gel strength can be adjusted by the time and intensity of the light;
- Add medium to the wells to cover the gel. Place the well plate in a 37°C incubator for 5 minutes. And then wash the sample and remove the medium;
- Add the cell suspension to the well plate. Change medium, observe, and photograph according to experimental design. (No special requirements for operation procedures).

Suggestions for 3D cell culture

- Cells were collected and resuspended in DexMA to prepare the cell suspension;
- Add cell suspension into the well plates;
(96-well plate: 50-100 μL/ well, 48-well plate: 100-300 μL/ well, 24-well plate: 300-500 μL/ well)
- Irradiate the wells with 405 nm light for 10-30 seconds to gelate, the gel strength can be adjusted by the time and intensity of the light;
- Add medium to the wells. Place the plate in a 37°C incubator for 5 minutes. And then wash the sample and remove the medium;
- Add fresh medium and incubate for a long time. Change medium, observe, and photograph according to experimental design. (No special requirements for

operation procedures).

Tips: Do not look directly at the light source.

