


Chitosan Methacryloyl (CSMA)

Product Component

Item	Concentration	Appearance	Solvent	Package Size	Remarks
CSMA solution	2% (w/v)	Light yellow clear solution	0.01M PBS Contains 0.25% (w/v) LAP	10 mL /bottle	Keep in dark

This instruction applies to EFL-S-CSMA-100K



CSMA molecular structure

Storage

room temperature, 1 months; 4°C, 3 months.

Period of validity

The date of manufacture is shown in the package.

Suggestions for cell culture

- The appropriate volume of CSMA solution was injected into the mold or cell culture well plate;
- 405nm wavelength ultraviolet lamp was used to irradiate CSMA for 10-30 seconds. Gel strength can be controlled by the duration and intensity of the light;
- The viscosity of the solution may rise during preservation, which is mainly caused by the non-covalent bond force. It can be restored to normal viscosity by high temperature autoclave procedure.

Tips: Do not look directly at the light source.

Scan the QR code on the right
for more information



WeChat Tencent